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PH-60/65

pH / mV / °C-meter

Inhoud / Content

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- 1 = BNC connector for pH sensor
- 2 = LCD display
- 3 = pH button
- 4 = Temperature button
- 5 = mV button
- 6 = Logo
- 7 = Temperature connector
- 8 = Temperature adjustment
- 9 = pH-4 adjustment
- 10 = pH-7 adjustment

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3.0

User manual:

- I To switch on the instrument press key marked "pH", after 3 min. the instrument will automatically shut off.
- II pH-key = 0.00-14.00 pH, mV = +/- 1000 mV, °C = -20/+120°C.
- III Set the temperature for the test sample by holding down the "°C" key while turning the "°C" dial on the right side to the correct temperature. The window automatically displays the pH value when the "°C" key is not pressed.
- IV When you have the pH-meter model PH-65, step III will be lapsed.
- V Soil, sand, peat and similar materials must be mixed with distilled water with a ratio of 2 parts water to 1 part soil. The C level on the sample vessel indicates 2 parts water and adding sample materials to the A level gives a ratio of 1 : 2. Water and other liquids may be measured directly.

Important for an accurate measurement:

- I Before measurement remove the electrode cover. This cover contains KCL preservation fluid which keeps the measuring part of the electrode constantly wet. This is very important for prolonging the live of the electrode and to prevent it for drying out.
- II Rinse the electrode gently in the water sample. The electrode is extremely fragile and care must be taken to avoid firm impact with stones and other hard materials.
- III Calibration must be carried out at least one time every 2 weeks, using a fresh calibration sample (pH-4 and pH-7). These can be obtained from your dealer or from Nieuwkoop. The "pH-4" and "pH-7" set screw on the left side of the instrument is used to calibrate the pH-meter. We advise to start your calibration always with the pH-7 liquid and then pH-4. Always rinse the electrode before you go from the pH-7 into the pH-4.
- IV After the measurement, rinse the electrode with distilled-water and store in the KCL-liquid.
- V When "LO BATT" appears in the display, change the battery, low battery power may cause failures.
- VI The pH-meter is also available as a combination instrument with EC and °C (option: automatic calibration) or as continuous measurement/controller.

3.1

Technical specifications:

Brand	:	Nieuwkoop Aalsmeer Holland
Measuring range	:	0-14 pH, +/- 1000 mV, -20,0/+120,0 °C
Ambient temperature	:	0-50 °C
Accuracy	:	0,05 pH
Resolution	:	0,01 pH
Electrode pH/°C	:	pH: epoxy gel, 100 cm cable with BNC-connector, °C: Pt1000
Power supply	:	1 x battery 9V 6F22
Measuring hours	:	± 250 hours
Dimension/weight	:	320x290x120 mm/ 2,1 kilos
Accessories	:	Carrier case, sample vessel, syringe, cleaning bottle, alcohol thermometer (model PH-60), distilled-water, calibration liquids (pH-4, pH-7), screwdriver and KCL storage liquid.
Warranty meter/electrode	:	1 year/ 3 months

The electrodes are filled with a gelled electrolyte solution and refilling is not necessary. When there is insufficient electrolyte solution, this type of electrode has probably been used at too high temperature and therefore replacement by another type of electrode is recommended.

Storage: When an electrode is to remain unused even for a small time, it is necessary to keep the glass-membrane in the KCL storage liquid in the storage reservoir.

Faults: Generally, faults are caused by:

- A) Sensitivity decrement of the membrane. When this occurs the electrode must be cleaned with a suitable detergent.
- 1) Depositions of hydroxides, lime iron hydroxide can be removed by immersing the electrode in a diluted solution of hydrochloric acid and then rinsing with water.
 - 2) Depositions of oil and fat can be removed by means of hot water together with domestic washing solution. When the result is unsatisfactory, a domestic abrasive may be carefully used.
 - 3) Albuminous depositions can be removed by means of a solvent containing: 8,5 ml concentrated hydrochloric acid, 10 gr. pepsin and 1 ltr. demi-water.
- Note: When polar solvents are used for special cleaning purposes, it is necessary to soak the electrode for some time after cleaning as the polar solvent influences the gel layer. When an a-polar solvent is used (benzene, ether, toluene) after treatment with a polar solvent (methanol, acetone) is necessary.
- B) Leakage to earth: Fault indications can be expected when the resistance between reference system and screening falls below 10/7 Ohm. When a fault occurs, first check if the electrode cable is in good condition and if the connector contacts and the terminals in the Nieuwkoop pH-meter or connection box are clean and dry.
- C) Short circuit by break or leakage: The pH-meter reading remains about pH-7 and is not dependant on the pH-value of the liquid being measured. If this occurs, the electrode should be replaced.
- D) Diaphragm poisoning: Strongly polluted liquids can cause the diaphragm to foul very quickly. As a result, there will be a transition resistance over the diaphragm. When this resistance exceeds 10/5 Ohm instability of reading and incorrect indication of the measured value will occur. Fouled electrodes can be cleaned with hot water together with a domestic washing solution if necessary. When fouling is caused by fat or hydroxides (lime) cleaning with an organic solvent or diluted acid respectively is recommended. It may be necessary with this non-flow electrode of the Nieuwkoop pH-meter, to place the electrode for a period in an electrolyte solution at 80 °C and let them remain until the solution has completely cooled.
- E) Poisoning: An electrode can be poisoned by the penetration of the process liquid or by diffusion of components of the process medium through the diaphragm. In addition, the inner liquid in combination with the dust from the process medium may give deposits which block the diaphragm (e.g. silver chloride + sulphide, deposits of silver sulphide. When poisoning occurs replacement of the electrode will normally be necessary. It is recommended if poisoning of a non-flow electrode occurs to use an alternate type of electrode.
- F) Poor conduction of the process liquid: When the electrical conduction of a process liquid is very low, an instability of reading will occur. This problem may be cured either by adding electrolyte solution to the process liquid or, if this is not possible, by making the reference input of the measuring instrument high impedance.
- G) Degrading of the glass membrane: Degrading has one of the following effects:
- Decrease of the speed of response
 - increase of the electrical resistance
 - decrease of the sensitivity
 - zero point shift.

The effect of ageing of the electrodes makes it necessary to carry out the re-calibration of the electrode system regularly. The frequency depends on the required accuracy, the composition of the sample to be measured and the temperature. E.g. if the lifetime of an electrode is 100% at the room temperature it will be 20% at 80 °C and only 5% at 120 °C. It is possible to reactive an aged electrode by immersing for 10 seconds in a PVC beaker containing a solution of vinegar (1 mol.) and potassium fluoride (1 mol.) Ratio 1:1. After this the electrode must be cleaned carefully.

Warning: As a result of a strong etching proportion of the acid, the handling should be done carefully and in consultation with a safety officer.

See further instructions of the pH-meter or combined pH/EC-meter.

Aalsmeerderweg 249s, 1432 CM
Postbus 78, 1430 AB
AALSMEER - HOLLAND
Tel. +31 (0)297 325836
Fax +31 (0)297 323167
E-mail: info@nieuwkoopbv.nl
Internet: www.nieuwkoopbv.nl

Uw leverancier/Your supplier:

--